

ES950



Born from the new generation of Orthoquad mobility scooters, the ES950 establishes a new bench mark for mid-size scooters. Powered by the new 950 watt 4 pole engine and coupled with two (2) 40 amp batteries, it delivers top performance with an impressive range. The ES950 has an amazing adjustable 4 wheel suspension for enhanced performance on rugged terrain.

It also has a long list of luxury standard features: full lighting package including directional signals, deluxe high back seat, enclosed front basket and delta tiller for easy throttle operation providing a safe mobility experience. Its compact design allows the ES950 to fit in all adapted public transportation. With the ES950, you get the best of everything. It is fast, powerful, comfortable and compact.

The ES950 will take you almost anywhere with style and safety.



Specifications

Length	48"	Charger	120V - 5A
Width	25"	Batteries	2 x 12V - 40Ah
Max. Speed	14.5 km/h	Wheels	10 "
Range (Up to)	30 - 40 km	Turning radius	1.65 m - 5,4'
Weight capacity	170 kg - 375 lbs	Braking distance	1.5 m - 5'
Suspension	Active front and rear	Drive Train	950W - 4 Pole
Transport weight W-O batteries	65 kg - 143 lbs	Max. climbing gradient	15 degrees