

# ORTHOQUAD

New Generation .COM

## RS1400

Orthoquad is proud to introduce its new generation full size scooter, the RS1400. Designed and developed for heavy duty performance it also includes all the luxury features that have come to be expected from the market. A long list of standard equipments as a full lighting package, deluxe high back seat and Delta tiller for easy throttle operation with fingers provided a safe and practical mobility experience. The RS1400 delivers a premium performance with 1400 watt four pole motor for exceptional power, torque and efficiency.

Excellent for outdoor performance with its 14" all-terrain tires, amazing ground clearance and a 4 wheel full suspension, the RS1400 will give you an optimal level of comfort on any surface. With the RS1400, you will find speed, comfort, range and reliability.

**The RS1400 is definitely your best choice!**



**RS1400**  
Dual seat edition

### Specifications

Length	63"	Charger	120V - 5A
Width	27"	Batteries	2 x 12V - 75Ah (group 24)
Max. Speed	19 km/h	Wheels	14 "
Range (Up to)	30 - 40 km	Turning radius	1.65 m - 5,4'
Weight capacity	180 kg - 400 lbs	Braking distance	1.8 m - 6'
Suspension	full active front and rear	Drive Train	1400W - 4 Pole
Transport weight W-O batteries	85 kg - 187 lbs	Max. climbing gradient	15 degrees